



Professional Tooth Whitening

Cosmetic tooth whitening

Application or first use only in the dental clinic

Method

Application time per session

Number of Applications



HOME WHITENING

Cosmetic tooth whitening treatment to be used at home by the patient with laboratory custom-made trays.

10%: overnight, min. 4 hours
16%: overnight or daytime, min. 2 hours

5 - 12
5 - 12



POWER WHITENING YF 6%

The fast method of cosmetic chairside tooth whitening in the dental clinic.

3 - 4 x 15 - 20 min

1 - 2 sessions



LIGHT WHITENING AC 6%

The most effective and fastest method of cosmetic chairside tooth whitening.

3-4 x 15 - 20 min

1 - 2 sessions

Medical tooth whitening

Application only with medical indication from dentist

Method

Application time per session

Number of Applications



IN-OFFICE WHITENING

Cosmetic tooth whitening treatment to be used at home by the patient with laboratory custom-made trays or as internal „walking bleach technique“ for endodontically treated teeth..

30 - 60 minutes in the waiting room

5 - 12



POWER WHITENING YF 40% (mixed 32%)

Medical tooth whitening chairside or for intracoronal whitening of endodontically treated teeth.

3 x 15 min (max. 4 x 15 min)

1 session, in exceptional cases 2 sessions



LIGHT WHITENING AC

The most effective and fastest method of medical chairside tooth whitening with light support. Also suitable for intracoronal whitening of endodontically treated teeth.

LW AC 22% HP
3 x 20 min (max. 4 x 20 min)

1 session, in exceptional cases 2 sessions

LW AC 32% HP
3 x 15-20 min (max. 4 x 15 min)

1 session, in exceptional cases 2 sessions